



# Ink Jet Printing of Organic Materials for Polymer Electronics

T. Wegener

GeSiM mbH, Bautzner Landstraße 45, D-01454 Großerkmannsdorf, Germany

E-Mail: wegener@gesim.de

## GeSiM Nano-Plotter

Pipetting device suitable for non contact pipetting and dispensing jobs in the sub micro liter range


### Key benefits:

- Excellent spot quality
- High reproducibility
- Low coefficient of variation from spot to spot
- Different spot sizes with one tip
- Arbitrary spot pattern with up to 16 tips
- Automatic function control of the pipettes
- Scalable workplate
- Head with 4 independent z-axis available


### Accessories:

- Humidifier
- Optical droplet detection
- Optical target recognition
- Coolable plate holder
- Coolable and removable slide tray
- Chiller
- Target visualisation
- Micro plate handler
- Z-level sensor

16 channel Nano-Plotter with video microscope dispensing on glass slides



16 channel Nano-Plotter with syringe pumps



## GeSiM micropipettes

Piezoelectric micro dispensers and pipettes with the ability to eject droplets down to 50 pico liters


Type	PicoTip	NanoTip	SIIP
Droplet diameter [µm]	40...70	80...100	100...125
Dispensing volume [pL]	33...180	250...520	600...1000

### Options:

- Different inlet configurations
- Integrated heater
- Layout adaption to specific dispensing tasks
- Multi nozzle print heads

### Successfully dispensed media:


- Organic solvents like toluene, DMF, dichlorobenzene, xylene, DMSO, ...
- Proteins (< 1 mg/ml), DNA
- Glycerol (50 % at rt, 100 % with heater)
- Ethylene glycol
- Up to 50 % PDMS
- Different adhesives
- Different liquid crystals



## Applications

### Targets and pattern recognition

#### Dispensing on an array of electrodes



The pictures show the finding of the spotting positions (left), the dispensing of the sample and the array with sample on the desired positions after spotting

The image recognition module of the NPC8/16 software contains a set of algorithms to detect spotting positions on several types of targets. For the electrode chip an algorithm which uses two rectangular alignment marks on the chip was chosen. These markers define a coordinate system and the user is free to set the spot positions anywhere by entering its coordinates with the help of our ChipDesigner program.

The pattern recognition algorithm finds the alignment marks and calculates the spot. Now the Nano Plotter™ is ready to spot exactly into the wells of the chip. After the camera returned to the chip we can observe the spotting result in a live video image. It is also possible to store images for documentation purposes. For better visibility just every second position is spotted in this example.

### Combinatorial Experimentation



Aspiration from a 96 micro well plate (left) and the stroboscopic images of tip 1 (Cy3, middle) and tip 2 (Cy5, right)

Fluorescence signals of Cy3 (left) and Cy5 (right)

One opportunity is to generate spots on a target consisting of defined amounts of different materials. To show this only two different materials are used here. The first material is deposited in horizontal lines with a constant number of droplets/spot within one line and an increasing number of droplets/spot from line to line. For the second material the number of droplets/spots is constant in the vertical lines. While travelling from one column to the other the number of droplets/spots has been increased by one. Following this procedure 400 spots with different composition have been fabricated.


The experiment has been performed using a two channel Nano-Plotter equipped with two Nano-Tips. The spotting position of the two tips has been adjusted using a video microscope and an image processing system.

The aspirated volume for both channels was 3 µl. A droplet volume of 320 pL has been used for both tips. The pitch of the resulting micro array was 600 µm. For the evaluation of the spots a fluorescence reader was used.

### Polymer Electronics

#### P3OT


- Line printing via overlapping of single dots
- Problems (flake like overlap) with trichloroethylene based solution due to too fast evaporation
- Well from chlorobenzene and xylene
- 150 .. 180 µm in width depending on concentration and jetting parameters with the given tip



oFET test devices with ink-jet printed line of dissolved P3OT: left: 1 mg/ml in trichloroethylen, right: 1 mg/ml in chlorobenzene

#### PEDOT:PSS

- test liquid: poly(styrenesulfonate)/poly(2,3-dihydrothieno[3,4-b]-1,4-dioxin), 2,8 wt% disp in water, elec. grade from Aldrich has been used
- the content of PEDOT and PSS in this solution is about 0,14% and 2,6 % respectively
- the viscosity is about 20 mPas.



Ink-jet printed lines and spots of PEDOT:PSS: left: lines dispensed at 1 kHz, width about 170 µm right: 10 x 10 array, spot size 145 µm

### High viscous liquids


#### Liquid crystals

##### Tested solutions:


- Merck Licristal ZLI-2222-100 (viscosity=12 mm<sup>2</sup>/s @ 23 °C)
- Merck Licristal MLC-6681 (viscosity=96 mm<sup>2</sup>/s @ 23 °C)

- Test solutions cover a wide viscosity range
- Although heating is recommended both liquids can be dispensed even at room temperature
- Droplet weight adjustable between 150 and 600 ng
- The maximum dispensing frequency is about 2 kHz depending on the temperature and the viscosity of the liquid

#### Droplet speed and weight of Merck Licristal ZLI-2222-100 (left) and MLC 6681 (right) both @ 43 °C



#### Ejected droplet of ZLI-2222-100



#### ZLI-2222-100 on glass pitch 350 µm

